

< Technic I >

Wooden Stick Preparation





Part facing diamond Adjust based on this part

Part that picks up the diamond


Diamond Nail Art I < Side Technic>


Tiara-Line

Items to Prepare


9 1 carat Nail Art Diamonds Thick soak-off gel Stainless steel mesh Aluminum cup

Items to prepare prior to completion of colored gel.





Looking at the Nail Bed from the front, the set face of the diamond should not be visible from the side of the diamond.



- Thick soak-off gel coating: tiny amount (when applied to nails, allow the UV gel to XXX degree)
- 2 Start from the center of the cuticle line.
- From the center of the cuticle line, at 1 mm gap away, place the 0.05ct diamond, and set its position.
- 4 After determining the 0.05ct diamond position, continue to place the remaining 8 diamonds (total of 9 diamonds), interchanging from left to right based on the sequence indicated on the left.
- After positioning all the diamonds according to the sequence, adjust the diamonds to match the curvature of the cuticle line. For confirmation, block the light with hands so that lighting condition is at a minimum.
- 6 allow the line to stay temporarily (shine light for 15 seconds)
- 7 The first layer Coating, coats the diamonds with the UV gel(BLUE)

I shine light for 15 seconds I

8 The second layer Coating, coats the entire surface circumference of the diamonds in the direction from 0.01 to 0.05 (GREEN)

, shine light for 15 seconds \downarrow

9 The third layer Coating, coats the entire surface circumference of the diamonds in the direction from 0.003 to 0.05(**ORANGE**)

shine light for 15 seconds I

10 The final Over Coating (PINK)

shine light for 30 seconds L

Use the acetone to clean off any UV gel that has not solidified.

IMPORTANT

Looking at the nail bed from the front with naked eyes, confirm that none of the diamonds are protruding from the sides, [especially these three points: the surface circumference of the diamond, sides and nail bed must all be confirmed].





